



CliniShield[®]
ADVANCE

360° Foam & Spray Skin Cleanser

CliniMed's cost-effective 'Clini' Range consists of CliniShield[®] Advance Barrier Film, CliniShield[®] Advance Protective Barrier Cream, CliniShield[®] Advance Protective Ointment, CliniShield[®] Advance 360° Foam & Spray Skin Cleanser and the CliniPeel[®] Medical Adhesive Removers.

CliniMed[®]

CliniShield[®] ADVANCE

360° Foam & Spray Skin Cleanser

Product Description

CliniShield Advance 360° Foam and Spray Skin Cleanser is a specially formulated non-rinse cleanser providing maximum protection from Moisture-Associated Skin Damage (MASD). It is a gentle, pH balanced cleanser which will minimise irritation and dryness of the skin, while moisturising.

Indications for use

CliniShield Advance 360° Foam and Spray Skin Cleanser indicated for use on intact and injured skin, removes urine and/or faeces from the skin following episodes of incontinence whilst minimising the risk of further skin damage. It effectively breaks down and removes stubborn dried bodily fluids like faeces and residue from skin barrier protectants. The product may be used on adults and children.

Benefits

- Cost effective: Cleanses and moisturises the skin
- The CliniShield Advance 360° Foam and Spray Skin Cleanser, can be used upside down, making it ideal for cleaning difficult-to-reach areas
- Can be used on intact or injured skin
- Non-rinse formulation minimising skin rubbing and potential friction damage
- Effectively removes dried-on bodily fluids
- Skin friendly, non-sting and non-sticky formulation
- pH balanced
- Spray & foam modes for easy application
- Alcohol, fragrance, parabens and phthalates free. Not made with natural rubber latex
- Even particle size distribution with non-drip dispenser
- Suitable for use on children

Contraindications

CliniShield Advance 360° Foam and Spray Skin Cleanser is not intended for direct use on mucosal membranes and broken skin. Do not use if known allergy or sensitivity. If in doubt, patch test a small area prior to use.

Caution

Use only as directed. For external use only. Do not use near and do not drip fluid into eyes, ears, nose, or mouth. For single patient use. May not be suitable for use on babies under 30 days old. Dispose of safely after use.

Instructions for use

1. To use, unlock by turning the nozzle to the spray position.
2. Firmly pull the trigger to dispense foam with the nozzle screen in closed position. Foam mode may be preferable for heavy soiling episodes and difficult-to-reach areas.
3. Gently open the nozzle screen to dispense spray. Spray mode might be preferable for lighter soiling episodes.
4. Apply directly to the skin or onto a washcloth. The 360° Foam and Spray Skin Cleanser can be used upside down, making it ideal for cleaning difficult-to-reach areas.
5. Gently wipe the area clean. No rinse required.
6. Pat dry or leave to air dry where possible.

CliniShield® ADVANCE

360° Foam & Spray Skin Cleanser

Warning



















Keep out of reach of children. Do not use in cases where there is a known sensitivity to any of the ingredients. Discontinue if irritation occurs and report to the manufacturer quoting the LOT number. Should redness or other signs of irritation appear, or irritation persists, consult a medical professional. If eye contact/ingestion occurs, seek immediate medical attention. Do not use if the packaging is damaged. In the event that the bottle is received damaged before use, contact CliniMed. Once open, products can be used through to expiry date. Keep out of reach of children. Avoid contact with eyes, if this happens, rinse thoroughly with water. Recycle empty bottles in accordance with local regulations. Trigger is not recyclable.

Ingredients

Aqua, Glycerine, PEG-6 Caprylic, Poloxamer 188, Coco-Glucoside, Phenoxyethanol, Chlorphenisin, Citric Acid

Storage and Shelf Life

Store away from direct sunlight at ambient temperature and humidity.
For shelf life, refer to the expiry date on the bottle.

	Manufacturer		Date of Manufacture/Country of Manufacture		Not Made with Natural Rubber Latex
	EU Representative				Recycling Label (PP)
	Catalogue Number		Use-by Date		Recycling Bottle (HDPE)
	Batch Code		Caution/Warning		Trigger Material Not Recyclable
	Medical Device		Keep Away from Sunlight		Do not use if packaging is damaged and consult instructions for use
	Unique Device Identifier		Keep Dry		
	Instructions for Use		Single Patient Multiple Use		

Product Code	Description
3832	CliniShield® ADVANCE 360° Foam & Spray Skin Cleanser 250ml

